


# MTT Cell Proliferation Kit

## INSTRUCTION MANUAL



**OZBIOSCIENCES**  
The art of delivery systems

[www.ozbiosciences.com](http://www.ozbiosciences.com)

# MTT – Cell Proliferation Kit

## Instruction Manual

Ready-to-use assay system for measurement of viability and proliferation rates based on spectrophotometric quantification of mitochondrial activity in living cells.

### **Catalog Number: MT01000**

For all other additional information, do not hesitate to contact our dedicated technical support ([tech@ozbiosciences.com](mailto:tech@ozbiosciences.com)).

#### **OZ Biosciences SAS**

163 avenue de Luminy  
Case 922, zone entreprise  
13288 Marseille cedex 09 - FRANCE  
Ph: +33 (0) 486 948 516  
Fax: +33 (0) 486 948 515  
[contact@ozbiosciences.com](mailto:contact@ozbiosciences.com)  
[order@ozbiosciences.com](mailto:order@ozbiosciences.com)

#### **OZ Biosciences INC**

4901 Morena Blvd,  
Suite 501  
San Diego CA 92117 - USA  
Ph : + 1-858-246-7840  
Fax : + 1-855-631-0626  
[contactUSA@ozbiosciences.com](mailto:contactUSA@ozbiosciences.com)  
[orderUSA@ozbiosciences.com](mailto:orderUSA@ozbiosciences.com)


[www.ozbiosciences.com](http://www.ozbiosciences.com)

# 1. Technology

## 1.1. Description

Congratulations on your purchase of the MTT – Cell Proliferation Kit!

The MTT – Cell Proliferation Kit is a colorimetric assay for measuring the mitochondrial reductase activity in living cells that reduce MTT to formazan dyes giving a blue/purple color. It is based on the cleavage of membrane-permeable yellow tetrazolium salt MTT to formazan crystals by metabolically active cells (figure 1).



**Figure 1:** MTT metabolization to formazan salt in metabolically active cells

A solubilization solution is then added to dissolve formazan into a colored solution. Spectrophotometric measurement of MTT-formazan at 570 nm allows quantitation of cell viability. Reagents used yield low background absorbance; a strong correlation between cell number and signal produced exists, allowing an accurate measurement of cell viability.

**This kit is designed for spectrophotometric quantification of cell growth, viability and proliferation and can be used as a direct indicator of cytotoxicity and apoptosis.**

This MTT – Cell Proliferation Kit is: Fast and Easy, Ready-to-use, Accurate, Economical and Stable under storage conditions.

## 1.2. Kit Contents

The kit contains sufficient reagents to perform 1000 assays in 96-well plate format.

Components	Quantity	Storage
MTT reagent 10X	10 x 1 mL	-20°C
Solubilization Solution (ready to use)	1 x 100 mL	4°C

## Stability and Storage

### Storage

Upon receipt and for long-term use, store all reagent tubes at the indicated storage conditions (see table above). Kit's components are stable for at least 1 year at the recommended storage temperature.

### Shipping condition

Room Temperature

## 2. Applications and Protocols

### 2.1. Usage

1. Plate the cells from 250 to 100 000 per well\* and follow your experimental protocol.
2. Perform the MTT assay to determine cell viability.
3. Incubate 30 minutes to 4 hours until purple crystals are formed
4. Solubilize precipitates and measure absorbance at 570 nm.
5. Determine viability.

\* number of cells depend on cell type (adherent/suspension), proliferation rate and metabolic activity.

### 2.2. Reagent preparation

- Dilute the 10X stock solution of MTT with sterile 1X PBS (pH 7.4). It is recommended to prepare fresh working solution each time the assay is performed.
- Excess of 1X solution can be returned at 4°C for short term storage (week) or -20°C for medium term storage (< 3 months).

Resulting solution should have a bright yellow color. Avoid repeated freeze/thaw cycles.

### 2.3. General Protocol for adherent cells in 96-well plate format

NOTE: 96 multiwell plates are the optimal vessel type for MTT method. Nevertheless it can be adapted to any vessel format, please refer to the table below.

1. Seed cells in a 96-well plate under standard culture conditions.
2. Carry out experiment by adding chemical compounds or biological agents to cells.

NOTE: at least two wells of cells should be kept untreated for viability control.

3. Wash cultured cells with 37°C pre-heated PBS.
4. Prepare **MTT working solution** as described in section 2.2.
5. Remove PBS and add 100 µL of **MTT working solution** to each well.
6. Incubate from 30 min to 4 h at 37°C (time depends on cell type/density/activity).

7. At the end of incubation time, add 100  $\mu$ L of **solubilization solution**. Pipette up and down several times to make sure the converted dye is completely dissolved.
8. Measure the absorbance of the converted dye on a plate reader at 570 nm and a reference wavelength at 650 nm.
9. Calculate the signal sample:  $OD_{570}-OD_{650}$  and express viability as a percentage of control cells.

#### 2.4. General Protocol for suspension cells

1. Seed cells in a 96-well plate under standard culture conditions.
2. Carry out experiment by adding chemical compounds or biological agents to cells.

NOTE: at least two wells of cells should be kept untreated for viability control.

3. Prepare **MTT working solution** as described in section 2.2.
4. Directly add 100  $\mu$ L of **MTT working solution** to each well.
5. Incubate from 30 min to 4 h at 37°C (time depends on cell type/density/activity).
6. At the end of incubation time, add 100  $\mu$ L of **solubilization solution**. Pipette up and down several times to make sure the converted dye is completely dissolved.
7. Measure the absorbance of the converted dye on a plate reader at 570 nm and a reference wavelength at 650 nm.
8. Calculate the signal sample:  $OD_{570}-OD_{650}$  and express viability as a percentage of control cells.

#### 2.5. General Protocol for larger samples

NOTE: This protocol is suitable to test viability on aliquots of trypsin-suspended adherent cells from large experiments or on suspension cells in large volumes.

1. Culture adherent or suspension cells as required.
2. Carry out experiment by adding chemical compounds or biological agents to cells.
3. For adherent cells: after trypsinization and washing, re suspend adherent cells. Take an aliquot of re-suspended adherent cells or suspension cells.
4. Transfer 250  $\mu$ L of each sample in a 1.5 mL tube
5. Prepare **MTT working solution** as described in section 2.2.
6. Directly add 250  $\mu$ L of **MTT working solution** to each tube.
7. Incubate from 30 min to 4 h at 37°C (time depends on cell type/density/activity).
8. At the end of incubation time, add 250  $\mu$ L of **solubilization solution**. Pipette up and down several times to make sure the converted dye is completely dissolved.
9. Centrifuge at 13.000 rpm for 2 min and transfer 100  $\mu$ L of supernatant in two wells of a 96-well plate.
10. Measure the absorbance of the converted dye on a plate reader at 570 nm and a reference wavelength at 650 nm.
11. Calculate the signal sample:  $OD_{570}-OD_{650}$  and express viability as a percentage of control cells.

**Volumes of solution recommended for various culture dishes are listed in the subsequent table.**

Type of culture dish	1X MTT working solution	Solubilization Buffer
96-well plate	100	100
24-well plate	250	250
12-well plate	500	500
6-well plate	1000	1000

### 3. Related Products

Description
<b>MAGNETOFECTION TECHNOLOGY</b>
Super Magnetic Plate ( <i>standard size for all cell culture support</i> ) Mega Magnetic plate ( <i>mega size to hold 4 culture dishes at one time</i> )
<b>Transfection reagents:</b>
PolyMag Neo ( <i>for all nucleic acids</i> )
Magnetofectamine™ ( <i>for all nucleic acids</i> )
NeuroMag ( <i>dedicated for neurons</i> )
SilenceMag ( <i>for siRNA application</i> )
<b>Transfection enhancer:</b>
CombiMag ( <i>to improve any transfection reagent efficiency</i> )
<b>Viral Transduction enhancers:</b>
ViroMag ( <i>to optimize viral transduction</i> )
ViroMag R/L ( <i>specific for Retrovirus and Lentivirus</i> )
AdenoMag ( <i>for Adenoviruses</i> )
<b>LIPOFECTION TECHNOLOGY (LIPID-BASED)</b>
Lullaby ( <i>siRNA transfection reagent</i> )
DreamFect Gold ( <i>Transfection reagent for all types of nucleic acids</i> )
VeroFect ( <i>for Vero cells</i> )
FlyFectin ( <i>for Insect cells</i> )
<b>i-MICST TECHNOLOGY</b>
Viro-MICST ( <i>to transduce directly on magnetic cell purification columns</i> )
<b>3D TRANSFECTION TECHNOLOGY</b>
3Dfect ( <i>for scaffolds culture</i> ) / 3DfectIN ( <i>for hydrogels culture</i> )
<b>RECOMBINANT PROTEIN PRODUCTION</b>
HYPE-5 Transfection Kit ( <i>for High Yield Protein Expression</i> )
<b>PROTEIN DELIVERY SYSTEMS</b>
Ab-DeliverIN ( <i>delivery reagent for antibodies</i> ) Pro-DeliverIN ( <i>delivery reagent for protein in vivo and in vitro</i> )
<b>PLASMIDS PVECTOZ</b>
pVectOZ-LacZ / pVectOZ-SEAP / pVectOZ-GFP / pVectOZ-Luciferase
<b>ASSAY KITS</b>
Bradford – Protein Assay Kit MTT cell proliferation kit β-Galactosidase assay kits (CPRG/ONPG)
<b>BIOCHEMICALS</b>
D-Luciferin, K <sup>+</sup> and Na <sup>+</sup> 1g X-Gal powder 1g / G-418, Sulfate 1g

Do not hesitate to contact us for all complementary information and remember to visit our website in order to stay inform on our last breakthrough technologies and updated on our complete product list. <http://www.ozbiosciences.com>

## Purchaser Notification

### Limited License

The purchase price paid for the MTT – Cell Proliferation kit by end users grants them a non-transferable, non-exclusive license to use the kit and/or its separate and included components (as listed in the section 1, Kit Contents). This reagent is intended **for internal research only** by the buyer. Such use is limited to the use in the product manual. Furthermore, research only use means that this kit and all of its contents are excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences.

Purchasers may terminate this License at any time by returning all MTT – Cell Proliferation Kit material and documentation to OZ Biosciences, or by destroying all MTT – Cell Proliferation Kit components. Purchasers are advised to contact OZ Biosciences with the notification that a MTT – Cell Proliferation Kit is being returned in order to obtain a refund and/or to expressly terminate a research only license granted through the purchase of the kit(s).

This document covers in full the terms of the MTT – Cell Proliferation Kit research only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

### Product Use Limitations

The MTT – Cell Proliferation Kit and all of its components are developed, designed, intended, and sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the handling of the kit components by following appropriate research lab practices.

For more information, or for any comments on the terms and conditions of this License, please contact:

Director of Business Development  
OZ Biosciences SAS  
Parc Scientifique et Technologique de Luminy  
Bâtiment Grand Luminy technopole  
Case 922, zone entreprises  
13288 Marseille Cedex 9, France  
Ph: +33 (0)4.86.94.85.16  
Fax: +33 (0)4.86.94.85.15  
E-mail: [contact@ozbiosciences.com](mailto:contact@ozbiosciences.com)

# CONTACTS

OZ Biosciences SAS  
163 avenue de Luminy  
Case 922, zone entreprise  
13288 Marseille cedex 09  
FRANCE

Ph: +33 (0) 486 948 516  
Fax: +33 (0) 486 948 515

[contact@ozbiosciences.com](mailto:contact@ozbiosciences.com)  
[order@ozbiosciences.com](mailto:order@ozbiosciences.com)  
[tech@ozbiosciences.com](mailto:tech@ozbiosciences.com)

OZ Biosciences INC  
4901 Morena Blvd,  
Suite 501  
San Diego CA 92117  
USA

Ph : + 1-858-246-7840  
Fax : + 1-855-631-0626

[contactUSA@ozbiosciences.com](mailto:contactUSA@ozbiosciences.com)  
[orderUSA@ozbiosciences.com](mailto:orderUSA@ozbiosciences.com)  
[techUSA@ozbiosciences.com](mailto:techUSA@ozbiosciences.com)

[www.ozbiosciences.com](http://www.ozbiosciences.com)

Follow us!

