


MULTIPLIX SHOULDN'T BE COMPLEX


Azure Imaging Systems

The Next Generation

Designed for easy Total Protein Normalization, Azure Imaging Systems bring even greater performance and streamlined user interface to the next evolution of quantitative Western blot imaging.




Simultaneous detection of total protein content with AzureRed and two-color NIR Western Blot using AzureSpectra secondary antibodies.




Digital image of 3-color Western Blot using Azure Biosystems 600 imager.


Chemi with
Azure Red for TPN



4 Channel Fluorescent
Western Blot Imaging



User Interface Designed
for Quantitative Westerns



Azure Imaging Systems

A complete, flexible imaging system line, designed to provide the flexibility needed for your applications.

Resolution and Applications	200	280	300	400	500	600
Resolution	5.4 MP	6 MP	9 MP	9 MP	9 MP	9 MP
UV Gel Imaging	✓	✓	✓	✓	✓	✓
Trans White Imaging	✓	✓	✓	✓	✓	✓
EPI Visible Color Imaging	✓	✓	✓	✓	✓	✓
Chemiluminescent Westerns		✓	✓	✓	✓	✓
Multiplex Visible Fluorescent Westerns				✓		✓
Multiplex NIR Fluorescent Westerns					✓	✓
Channel for TPN			Optional	✓	Optional	✓

Aurogene S.r.l.

Some of our customers



UNIVERSITÀ
DEGLI STUDI
DI MILANO



Fondazione I.R.C.C.S.
Istituto Neurologico Carlo Besta



Auxologico
Ricerca e cura per la tua salute
IRCCS



UNIVERSITÀ
DEGLI STUDI
DI TRIESTE



UNIVERSITÀ
di VERONA



evotec



ALMA MATER STUDIORUM
UNIVERSITÀ DI BOLOGNA



UNIVERSITÀ
DEGLI STUDI
DI TERAMO



SAPIENZA
UNIVERSITÀ DI ROMA



SISTEMA SANITARIO REGIONALE
IRCCS
LAZZARO SPALLANZANI



TOR VERGATA
UNIVERSITÀ DEGLI STUDI DI ROMA



Aurogene S.r.l.

Via dei Lucani 51,53,55

Phone: 06 98185510

Mail: info@aurogene.eu

www.aurogene.eu

Aurogene S.r.l.